

EXPERIMENT 1

Gamma-Ray Absorption and Counting Statistics

Note: Please read the Radiation Safety Regulations at the back of this book

Introduction

The detecting medium in a scintillation counter is a solid or liquid material called a phosphor, which emits tiny flashes of light (scintillations) when struck by ionising radiations. The flashes are observed by a photomultiplier which produces an electrical pulse when a light flash falls on it.

In this experiment the phosphor is a 50.8 mm long \times 50.8 mm diameter cylindrical crystal of thallium-activated sodium iodide i.e. NaI(Tl). This material has a density of $3.67 \times 10^3 \text{ kg m}^{-3}$. Because of this and the high atomic number ($Z = 53$) of the iodine atoms, γ -rays interact strongly with the material, resulting in a γ -ray sensitivity very much greater than that of the Geiger-Müller counter. The crystal is enclosed in an airtight capsule of thin aluminium with diffusely-reflecting sides and a transparent base through which the light of the scintillations emerges. The base is in contact with the face of an 11-stage venetian-blind type photomultiplier.

It is arranged that no light except that from the scintillations can reach the photomultiplier and the crystal-photomultiplier combination is also surrounded by a lead shield which excludes unwanted background radiation (see Fig. 1.1).

The amount of light emitted in a scintillation is accurately proportional to the energy-loss of the particle producing the flash. The particles responsible for the flashes are the secondary electrons produced by the photons, not the photons themselves. When the phosphor is struck by γ -rays, for example, interaction takes place by both photoelectric and Compton effects. In the first case the photon loses all its energy to a photoelectron, which in turn yields its entire energy to the phosphor. The flash intensity in this case is proportional to the energy of the incident photon. In a Compton collision the photon transfers only a portion of its energy to the secondary electron; the secondary electron is brought to rest in the phosphor and causes the emission of a flash whose intensity is proportional to the energy-loss of the photon. Since in every case the photomultiplier produces pulses whose heights are proportional to the flash intensities, the pulse-spectrum resulting from the impact of monoenergetic γ -rays on the crystal is a mixture of pulses derived from full and partial energy-losses (see Figs. 1.2 and 1.7).

The pulses corresponding to full energy-loss, produced as they are by the photoelectric effect, fall




Figure 1.1: Schematic of photomultiplier and shielding




Figure 1.2: Pulse spectrum resulting from the impact of monoenergetic γ -rays on the crystal




Figure 1.3: Gamma-ray detection system

in a region called the photopeak while the pulses in the Compton ‘tail’ or continuum result from Compton collisions in which the scattered photon retains a portion of its energy. It is seen that even in the case of monoenergetic radiation the pulse-spectrum is not simple; when the incident radiation is heteroenergetic the pulse-spectrum may be very complex indeed.

More details of the processes of photoelectric absorption, Compton scattering and pair production may be found in the Appendix.

The Gamma-Ray Detection System

The radioactive source used is ^{137}Cs with a half life of 30 years. It was produced as a fission fragment by the neutron bombardment of ^{235}U . This isotope emits beta rays of energies 512 keV (94.6%) and 1.174 MeV (5.4%) that are absorbed in plastic before reaching the detector as well as the monoenergetic gamma rays of energy 662 keV.

The apparatus is set up as shown in the block diagram, Fig. 1.3. The high voltage power supply, preamplifier, main amplifier and analog to digital detector are all located in one box beside the computer and printer.

(CAUTION! Do not switch on until the demonstrator has checked the system).

The amplified photomultiplier pulses are digitised using an analog to digital converter and the digital pulse height is stored in computer memory. The computer memory is configured to place a count in any one of 1020 separate channels according to the digital value of the pulse height. The smallest pulse amplitude E is preset and the photopeak should occur somewhere between channels 500 and 700. The width of each channel fixes the amplitude range ΔE within which a pulse must fall if it is to be counted in that channel. The width of each channel is approximately 1 keV. When the photomultiplier produces a pulse whose amplitude is between E and $(E + \Delta E)$, the pulse height analyser increments by one the count in that channel. Using this method the full spectrum of pulse heights and hence gamma ray energies are recorded and stored. The width of the photopeaks is approximately 80 channels. The width is determined by the energy resolution ($\approx 12\%$) of the phosphor and not the monoenergetic gamma rays of energy 662 keV.

Instructions on the Operation of the Computer Programs and Computer Graphics

Remove any lead or iron absorbers from the gamma ray beam. Insert floppy disc containing the programs into the drive. Use the arrow keys to select (Note CR = carriage return)

- a) CONFIGURATION : FULL and a CR
- b) SAVE : ON and a CR
- c) HPRINT : ON and a CR

The graphics display has now been configured for the monitor and printer. Load the pulse height analysis (PHA) program after the computer has replied A> by typing PHA and a CR.

A number of options are now available

- a) Integration time. You may select a time between 5 and 180 seconds but normally 30 is sufficient. Type 30 followed by CR. The photopeak should occur between channel numbers 500 and 700 and the spectrum shown in Fig. 1.2 should be displayed on the screen.
- b) To adjust the position of the left marker on the screen use the keys ← and →. The channel number of the marker is displayed on top of the screen along with the count in that channel.
- c) To adjust the position of the right marker on the screen use the keys ↑ and ↓. The channel number of the marker is displayed on top of the screen along with the count in that channel.
- d) The total count between the two markers is displayed on the bottom of the screen containing the plot along with the separation between the two markers.
- e) The key F9 interrupts the program and you need to type CR for the program to continue.
- f) The key F1 resets the counters and starts the integration.
- g) You can obtain a copy of the spectrum on the screen by typing F9 followed by the shift key (⇧) and Prt SC key (i.e. upper case selection) and a CR.
- h) To input a new integration time, press the F5 key followed by a CR. When prompted you should type PHA followed by CR and input the new integration time.

You are now ready to continue the experiment.

You should obtain a printout of the spectrum on the screen using the procedure described in (g). Locate the position of the photopeak using either the left or right marker.

As shown in Fig. 1.4 the spectrum you have recorded is a differential spectrum, i.e. $N(E, E + \Delta E)$ versus E or the number of counts between E and $E + \Delta E$ as a function of the energy E .

You should readily convince yourself, by plotting the differential spectrum in the integral form i.e. $N(>E)$ versus E , that you cannot easily locate the photopeak. This can be done by adjusting the spacing between the cursors. The integral spectrum obtained from the differential spectrum, is also shown in Fig. 1.4.

You should also record the gamma ray spectrum of ^{60}Co and compare it with that for ^{137}Cs (see Figs. 1.7 and 1.8). The backscatter peak depends on the geometry of the shielding and may not be as prominent in your data as in Figs. 1.7 and 1.8.




Figure 1.4: Differential and integral spectra.

Experiment 1 - Linear absorption coefficients of Lead and Iron

The apparatus has been set up (see Fig. 1.3) and checked. The ^{137}Cs photopeak is now located by the method just described. Place the left hand marker at the position A in Fig. 1.2 and the right hand marker at position B. You can now determine the count in the photopeak and scattered photons are eliminated. Lead absorbers are now placed over the source and the count rate between A and B is recorded as a function of the absorber thickness.

When gamma-radiation passes through matter, it undergoes absorption primarily by Compton, photoelectric and pair production interactions. The intensity of the radiation is thus decreased as a function of distance in the absorbing medium. The mathematical expression for the intensity I is given by

$$I = I_0 e^{-\mu x}$$

where

I_0 = original intensity of the beam (no lead)

I = intensity transmission through the absorber to a distance, depth or thickness x

μ = linear absorption coefficient (cm^{-1}).

Plot the natural log of the counts as a function of absorber thickness. From the slope of the straight line determine μ (in units of cm^{-1}) and compute μ/ρ (in units of cm^2/g) where ρ is the density of lead. Compare your result with the value of μ predicted by the accompanying graph (Fig. 1.6) of the absorption coefficient versus gamma ray energy.

Repeat the above procedure for the iron absorbers and determine μ for iron.

Experiment 2 - Counting Statistics

Purpose As is well known, each measurement made for a radioactive sample is independent of all previous measurements, because radioactive decay is a random process. However, for a large number of individual measurements the deviation of the individual count rates from what might be termed the “average count rate” behaves in a predictable manner. Small deviations from

Experiment 1. Gamma-Ray Absorption and Counting Statistics

the average are much more likely than large deviations. In this experiment we will see that the frequency of occurrence of a particular deviation from this average within a given size interval can be determined with a certain degree of confidence. Fifty independent measurements will be made, and some rather simple statistical treatments of the data will be performed.

The average count rate for N independent measurements is given by

$$\bar{R} = \frac{R_1 + R_2 + R_3 + \dots + R_N}{N} \quad (1.1)$$

where R_1 = the count rate for the first measurement, etc., and N = the number of measurements

In summation, notation \bar{R} would take the form

$$\bar{R} = \frac{\sum_{i=1}^{i=N} R_i}{N} \quad (1.2)$$

The deviation of an individual count from the mean is $(R - \bar{R})$. From the definition of \bar{R} it is clear that

$$\sum_{i=1}^{i=N} (R_i - \bar{R}) = 0 \quad (1.3)$$

The standard deviation $\sigma = \sqrt{\bar{R}}$.

Procedure

1. Insert the 15.8 grm cm^{-2} lead absorber in the beam and adjust the left-hand and right-hand markers so that approximately 1000 counts can be recorded between the markers in a time period of 30 seconds.
2. Take 50 independent runs of duration 30 seconds and record the values in Table 1.1. (Note that you will have to extend Table 1.1; we have shown only ten entries). The count values R may be recorded directly in the table since for this experiment R is defined as the number of counts recorded in a 30 second interval.
3. With a calculator determine \bar{R} from Eq. 1.1. Fill in the values of $R - \bar{R}$ in Table 1.1. It should be noted that these values can be either positive or negative. You should indicate the sign in the data entered in the table.

Data Analysis:

Exercise a. Calculate σ , and fill in the values for σ and $(R - \bar{R})/\sigma$ in the table, using only two decimal places. Round off the values for $(R - \bar{R})/\sigma$ to the nearest 0.5 and record these values in the table. Note that in our Table 1 we have shown some typical values of $(R - \bar{R})/\sigma$ and the rounded-off values.

Exercise b. Make a plot of the frequency of the rounded-off events $(R - \bar{R})/\sigma$ vs the rounded-off values. Fig. 1.5 shows this plot for the ideal case of a normal distribution.

Table 1.1: Typical Values of $(R - \bar{R})/\sigma$ and $(R - \bar{R})/\sigma$ rounded off; listed for illustrative purposes only

Run	R	R - \bar{R}	$(R - \bar{R})/\sigma$		$(R - \bar{R})/\sigma$ (Round'd Off)	
			Typical	Measured	Typical	Measured
1			-0.15		0	
2			+1.06		+1.0	
3			+0.07		0	
4			-1.61		-1.5	
5			-1.21		-1.0	
6			+1.70		+1.5	
7			-0.03		0	
8			-1.17		-1.0	
9			-1.67		-1.5	
10			+0.19		0	

Note that at zero there eight events, etc. This means that in our complete rounded-off data in Table 1.1 there were eight zeros. Likewise, there were seven values of +0.5, etc. Does your plot follow a normal distribution similar to that in Fig. 1.5? If time permits, increase the amount of lead absorber and cursor positions to give a mean rate of $\approx 100s^{-1}$ and take 50 independent 30 second runs. Compute \bar{R} and σ . Compare your results for σ for the two different mean rates.

Question

1. Is it possible for a γ -ray to interact with a free electron by means of the photoelectric effect?

Appendix

Interactions of γ -Ray Photons with Matter

Fig. 1.9 shows some of the processes that can occur when a γ -ray photon enters a solid detector. The photon can Compton scatter several times; after each scattering, the photon loses some energy and a free electron is produced. Gradually the photon suffers either of two fates: it continues the repeated Compton scattering, eventually becoming so low in energy that photoelectric absorption occurs and the photon vanishes, or else it wanders too close to the edge of the crystal and scatters out of the crystal. The energy of the photon is converted into electrons (photoelectrons or Compton-scattered electrons), which have a very short range in the crystal, and which therefore lose energy very rapidly, by creating light photons in a scintillator or electron-hole pairs in a semiconductor detector. We can assume that all of this energy is absorbed, and we will refer to this quantity as the energy deposited in the detector by the original photon. If the original photon eventually suffers photoelectric absorption, the energy deposited is equal to the original γ -ray energy. If it scatters out of the crystal, the energy deposited is less than the original photon energy.




Figure 1.5: Typical Plot of Frequency of Rounded-Off Events vs the Rounded-Off Values




Figure 1.6: Narrow Beam Absorption Coefficients

Experiment 1. Gamma-Ray Absorption and Counting Statistics




Figure 1.7: NaI(Tl) Spectrum for ^{137}Cs




Figure 1.8: NaI(Tl) Spectrum for ^{60}Co

Experiment 1. Gamma-Ray Absorption and Counting Statistics




Figure 1.9: Processes occurring in γ -ray detection (1) The photon Compton scatters a few times and eventually leaves the detector before depositing all its energy. (2) Multiple Compton scattering is followed by photoelectric absorption, and complete energy deposition occurs. (3) Pair production followed by positron annihilation, Compton scattering, and photoelectric absorption; again, complete energy loss occurs. (4) One of the annihilation photons leaves the detector, and the γ -ray deposits its full energy less 511 keV. (5) Both annihilation photons leave the detector, resulting in energy deposition of the full energy less 1022 keV. Processes (4) and (5) occur only if the γ -ray energy exceeds 1022 keV.

Let's consider how much energy is given to the scattered electron in a single Compton event. The electron kinetic energy is given by:

$$T_e = E_\gamma - E'_\gamma = \frac{E_\gamma^2(1 - \cos \theta)}{mc^2 + E_\gamma(1 - \cos \theta)}$$

(see Krane, Introductory Nuclear Physics for the derivation).

E_γ is the energy of the incident γ -ray photon.

E'_γ is the energy of the Compton-scattered γ -ray photon.

Since all scattering angles can occur in the detector, the scattered electron ranges in energy from 0 for $\theta = 0^\circ$ to $2E_\gamma^2/(mc^2 + 2E_\gamma)$ for $\theta = 180^\circ$. These electrons will normally be totally absorbed in the detector, and (if the scattered photons escape) they contribute to the energy response of the detector a continuum called the *Compton continuum*, ranging from zero up to a maximum known as the *Compton edge*. The peak at $E = E_\gamma$ corresponding to complete photoelectric absorption (called the *full-energy peak or photopeak*) and the Compton continuum are shown in Fig. 1.10.




Figure 1.10: A typical response of a detector to monoenergetic γ -rays. The photopeak results from the γ -ray losing all its energy in the detector, as in events 2 and 3 in Fig. 1.9. The Compton continuum consists of many events of type 1 while the single- and double-escape peaks result from processes 4 and 5. The detector energy resolution might tend to broaden all peaks more than they are shown here, and multiple Compton scattering will fill in the gap between the Compton edge and the photopeak. The escape peaks appear only if the γ -ray energy is above 1.022 MeV.

We have so far neglected the third process of γ -ray interactions in the detector, that of pair production. The positron and electron are created with a total kinetic energy of $E_\gamma - 2mc^2$, and loss of this energy in the detector would result in a peak at the full energy. However, once the positron slows down to an energy near that of the atomic electron, *annihilation* takes place, in which the positron and an atomic electron disappear and are replaced by two photons of energy mc^2 or 511 keV. These two photons can travel out of the detector with no interactions, or can be totally or partially absorbed, through Compton scattering processes. We therefore expect to see peaks at $E_\gamma - 2mc^2$ (when both photons escape), $E_\gamma - mc^2$ (when one escapes and the other is totally absorbed) and E_γ (when both are totally absorbed). These single- and double-escape peaks are shown in Fig. 1.10.

The relative amplitudes of the photopeak, Compton continuum, and escape peaks depend on the size and shape of the detector. In general, the larger the detector, the smaller is the Compton continuum relative to the photopeak, for there is a smaller chance of a Compton-scattered photon surviving from the center to the surface without interacting again. Similarly, in a large detector there is a greater chance of capturing one or both of the 511-keV annihilation photons.

Taken from "Introductory Nuclear Physics", Kenneth R. Krane.

Reference

1. 'Introductory Nuclear Physics', Kenneth R. Krane

WWW :

- Radiation Physics and Radiocarbon Research Lab in UCD - <http://www.ucd.ie/radphys>
- Gamma-ray Astronomy - <http://imagine.gsfc.nasa.gov/docs/introduction/gamma-information.html>
- European Gamma-Ray Mission 'INTEGRAL' - <http://astro.estec.esa.nl/Integral>